







ID	Title	External_doc	External_risks	External_validation	External_validation_test
 100208-5236	Anbindung 3G/4G Modem über Ethernetschnittstelle	<p>Base</p> <p>Internal GSM modem (2G) might be replaced by external gateway with 3G/4G connection.</p> <p>2G modem will be deactivated if external modem is selected.</p> <p>External modem supports SMS alarming in 3G/4G communication, but with restricted functionality.</p>	<p>- SMS size generally limited to max. 140 char (for internal modem too)</p> <p>- no supervision of SMS delivery possible</p>	yes	Alarm functionality by SMS has to be checked if external modem is used instead of internal GSM module.
 100208-5337	FW - Base hängt beim Herunterfahren "Shutdown - Loading CommUp"	<p>Base V2</p> <p>in some cases the Saveris Base was blocked during shutdown.</p> <p>The issue has been solved.</p>	none		
 100208-5061	Alarm - AM mismatch between colors and values in table and alarmlog	<p>Ethernetprobe</p> <p>In some cases a mismatch between colours and values in table and alarmlog was observed.</p> <p>Also the behaviour between alarm by email and alarm by SMS was different.</p> <p>Now the behaviour is synchronized for email and SMS alarms - the corresponding colour fits to the listed start and end of alarm timepoints in both cases.</p>	none		
 100208-5169	Unnötiger Eintrag in der Konfigurationsvorgabe der Ethernet-Logger	<p>Ethernetprobes</p> <p>In the config memory range of ethernetprobes a value was stored which is only in use for Radio</p>	none		

		<p>probes. This config-value for radio probes will not be stored to ethernetprobes furtheron.</p>	
 100208-5112	Alarm - Dauer Alarmverzögerung wird von Türkontakt-Messpunkten verkürzt	<p>Radio probes V1 and V2 Alarm-delay is based on numbers of standard measuring cycles . If alarm delay is used in combination with door-contact, the alarm comes earlier than expected because door-opening creates additional timestamps between the standard cycle timestamps and reduces the alarm-delay. Now this reduction is blocked.</p>	none
 100208-4749	Alarm - Neue Darstellung der Alarmnachrichten	<p>The e-mail alarm message body has been re-formatted and enhanced. Here is a sample:</p> <pre> TestoSaveris-Alarm 67373337-1 - //nnpce1610/tsaveris ----- AlarmID: 67373337 Condition: 18.2 °C < 20.0 °C Date Time: 19.09.2016 13:21:00 -----Source----- <devicename> SN: <serialnumber> Channel: <channelnumber> Zone(s): <Zone>[,morezones] Project: <projectname> </pre>	none

<Alarm message text>
