

Temperature probe – conformity-evaluated

For use in pathology,
forensic medicine and laboratories.

Spherical measuring tip

Probe shaft and connecting cable can be customized

Compatible with the conformity-evaluated
temperature measuring instrument testo 112



This conformity-evaluated temperature probe, with spherical measuring tip to protect the tissue when inserting, is ideal for temperature measurements in pathology, forensic medicine or laboratories. The probe is compatible with the conformity-evaluated temperature measuring instrument testo 112.

The probe shaft is available in various lengths and diameters. The length and material of the connecting cable can also be customized. The probe is also available as a Pt100 variant.

Temperature probe - conformity-evaluated

Temperature probe similar to 0614 1212 01, conformity-evaluated, but with probe shaft \varnothing 4 mm, length 150 mm. Probe shaft made of stainless steel 1.4571 with spherical measurement tip. Plastic handle (0400 9736 weis). fixed connection cable 1.6 m with PUR insulation, plug fits testo 112.



Order no. 0600 9999 01
ID no. 0699 4471/4

Temperature measuring instrument testo 112

testo 112, 1-channel temperature measuring instrument NTC/Pt100, conformity-evaluated, includes battery



Order no. 0560 1128 01

Standard version

Sensor type	NTC
Measuring range	-25 to +120 °C
Accuracy	± 0.5 % of m.v. (+100 to +120 °C) ± 0.2 °C (-25 to +74.9 °C) ± 0.4 °C (remaining meas. range)
Handle material	Plastic
Probe tube material	Stainless steel 1,4571
Measuring tip	Spherical
Probe tube length	150 mm
\varnothing probe tube	\varnothing 4 mm
Connector	8-pin mini DIN connector
Connecting pipe	Length 1.6 m, PUR-insulated

Possible variants

Handle material	Stainless steel
Sensor type	Pt100 class A
Probe tube length	100 to 300 mm
\varnothing probe tube	3 to 5 mm
Cable length	Up to max. 5 m
Cable insulation	FEP

We are happy to adapt probes to suit your individual requirements.
Simply send us a short mail to: vertrieb@testo.de.



Testo SE & Co. KGaA
Celsiusstrasse 2, D-79822 Titisee-Neustadt, Germany
Phone +49 7653 681-700
Fax +49 7653 681-701
vertrieb@testo.de